



Recipients: Veterinarians from the swine sector
Sender: Diagnostic Service
Date: 2018-11-29
Objet: **Enzyme change for PRRSV detection test by PCR**

Dear customers,

We would like to inform you that the Molecular Diagnostic laboratory of the Diagnostic Service has recently changed the reagents used in the PRRSV testing by PCR. This change allows an increase of sensitivity of the test. The detection limit thus changed from 44 genome copies to about 3 genome copies.

Quantification was adjusted for this change. You could therefore observe a positive result with a lower quantification of about 1 LOG for either tissues, sera, oral fluids or other liquids.

We are confident that this change will meet your needs.

If you have any questions, do not hesitate to contact us.

Véronique Allard
Technical and Customer Service Manager

Diagnostic Service
Faculty of Veterinary medicine
Université de Montréal
veronique.allard@umontreal.ca
450-773-8521 ext. 8670

Customer Service

Phone: 450-778-8151 (direct line)
450-773-8521 ext. 8243

Fax: 450-778-8107

service-diagnostic@medvet.umontreal.ca
www.servicedediagnostic.com